

About Signia

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath. With Signia, you can Be Brilliant.

Find out more on signia-pro.com

Or ask your hearing care professional for more information on Signia hearing aids, useful accessories and the Signia app.



Compatible with Android and iOS. The app is free of charge.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by WSAUD A/S is under license. Other trademarks and trade names are those of their respective owners.

The 'Qi' symbol is a trademark of the Wireless Power Consortium.

Find us on

-  signia-pro.com
-  facebook.com/signiahearing
-  youtube.com/signiahearing
-  instagram.com/signiahearing

Order No. 0xxx-99C1-7600 | Printed in Germany
© 09.2023 WSAUD A/S

Be Brilliant™

Product Portfolio 2023 – 2nd Half





Be Brilliant™

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath.

With Signia, you can Be Brilliant.

Contents



Signia IX	4
Integrated Xperience	6
Silk Charge&Go IX	12
Pure Charge&Go IX / Pure Charge&Go T IX	13
Signia AX	14
Augmented Xperience	16
Styletto AX	18
Pure Charge&Go AX / Pure Charge&Go T AX	19
Insio Charge&Go AX	20
Pure 312 AX	21
Signia Xperience	22
Active / Active Pro	24
Motion Charge&Go X, P X, SP X	25
Styletto X	28
Silk X	29
Pure Charge&Go X	30
Pure 312 X	31
Signia Nx	32
Essential	36
Basic	42
Accessories	48
Hearing care professional support	56
App	58
Color overview	60
Feature overview	62

Signia IX

Unleash the power of conversation

Performance levels

- 7IX
- 5IX
- 3IX

Meet the world's first hearing aid with a multi-stream architecture capable of accurately pinpointing multiple moving speakers in real-time conversation.

Improved hearing is about more than just sound amplification. It is about empowering individuals to immerse themselves in their surroundings, create meaningful connections, and cherish every moment of life's experiences.

Leveraging our profound understanding of dynamic conversations, including the head movements of the wearer and changing positions of speakers, we crafted Signia Integrated Xperience. This groundbreaking technology ensures uninterrupted conversation engagement in noise, addressing a critical challenge many hearing aid wearers face.

Our commitment is to support wearers in leading active lives, unshackled from the limitations of traditional hearing aids and free to contribute and foster deeper connections.

Experience the transformative power of real-time, dynamic conversation enhancement with Signia Integrated Xperience.



Signia Integrated Xperience is available from 7 to 3 performance levels for Pure Charge&Go IX, Pure Charge&Go T IX and Silk Charge&Go IX. Now you can offer all your customers the transformative power of real-time, dynamic conversation enhancement with Signia Integrated Xperience.

ITE

New



Silk Charge&Go IX

World's First
Rechargeable
Instant fit CIC

RIC

New



Pure Charge&Go IX

New



Pure Charge&Go T IX





RealTime Conversation Enhancement

The science of keeping up with conversations

Signia Integrated Xperience's revolutionary RealTime Conversation Enhancement analyzes, augments, and adapts to the dynamic flow of conversations. This lets wearers enjoy a seamless and enriched hearing experience.



Analyze

Processes 192,000 data points every second to identify speakers and analyze conversation dynamics, authentically reflecting real-life conversations.



Augment

Uses data gathered to control multiple focus streams that lock on to and dynamically highlight each speaker's voice amidst the noise. Meanwhile, Augmented Focus technology[2] processes background sound to ensure a comfortable, immersive conversational experience without distractions.



Adapt

Our patented multi-stream architecture updates 1000 times per second, merging multiple focus streams into a dynamic auditory space that adapts to evolving conversations. This allows for real-time tracking and enhancement of multiple moving speakers, ensuring wearers can fully engage in group discussions.





Integrated Xperience feature overview

Signia Integrated Xperience, featuring RealTime Conversation Enhancement, is the heart of our innovative hearing aids. This technology empowers wearers to engage dynamically in conversations, ensuring every word is understood, every laugh appreciated, and every story fully experienced.

Features	Form Factor		Description	Performance levels		
	Pure	Silk		7 IX	5 IX	3 IX
Conversation management						
RealTime Conversation Enhancement	✓		Provides uninterrupted engagement in noisy group conversations by seamlessly tracking the locations of all conversation partners and enhancing their speech.*	●●●●●	●●●●○	●●○○○
Conversation Source Analyzer	✓		Constantly scouts the acoustic environment to identify number of talkers and their location.	✓	✓	✓
Conversation Dynamics Analyzer	✓		Analyzes conversation dynamics including turn-taking and timing, locations of pure noise without speech, and levels of background noise.	✓		
Focus Streams	✓		Isolate and track speech independently.	3	1	
Conversation Booster	✓		Dynamically enhances individual speakers.	●●●●●	●●●●○	●●○○○
Augmented Focus	✓		Separates the focus sound (speech) from surrounding sounds, creating a strong contrast.*	✓	✓	✓
Focus Stream	✓		Enhances speech signals from the front hemisphere.	1	1	1
Surrounding Stream	✓		Processes environmental sounds to preserve environmental awareness without masking speech.	1	1	1
Own Voice Processing 2.0 (OVP™ 2.0)	✓		Uses real-time recognition of the wearer's voice to deliver a natural own-voice impression.*	✓	✓	✓
Dynamic Soundscape Processing 3.0	✓	✓	Integrates inputs for precise hearing aid optimization, including (when available): acoustic sensors, motion sensors, OVP 2.0, noise detection, SNR analysis, and RealTime Conversation Enhancement.	●●●●●	●●●●○	●●●○○
SpeechStabilizer	✓	✓	Preserves the natural envelope of speech in any situation, including very noisy situations, which preserves acoustic cues that are important for speech understanding.	✓	✓	✓
Motion Sensor	✓		Detects whether the wearer is walking, running, or sitting still to guide how sound is optimized.	✓	✓	✓
e2e wireless 4.0	✓	✓	Controls coupling, synchronization, CROS/BiCROS, and a range of sound features.*	✓	✓	✓
Speech and Noise Management	✓	✓	Reduces listening effort using sophisticated sound processing.	✓	✓	✓
TruEar 360	✓		Reproduces benefits of the natural pinna effect by highlighting the difference between front and back sources of sound.	✓	✓	✓
Binaural Directionality	✓		Provides scientifically proven enhanced speech understanding in noise.*	✓	✓	✓
Binaural OneMic Directionality		✓	Microphones from each hearing aid work together to create a directional listening pattern preferentially amplifying speech from the front to support hearing in noisy environments.*	✓	✓	✓
Spatial SpeechFocus	✓		Steers directional microphone beam to the dominant speech source (left, right, or behind the wearer automatically in a car situation or in a dedicated "Strolling" program).*	●●●●●	●●●●○	

*) Available for bilateral fittings with e2e wireless link
**) Requires Bluetooth connection, not available in Silk Charge&Go IX

Features	Form Factor		Description	Performance levels		
	Pure	Silk		7 IX	5 IX	3 IX
Listening Comfort						
Feedback Cancellation	✓	✓	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓
eWindScreen™	✓	✓	Reduces the interference of wind noise when outdoors.	●●●●●	●●●●●	●●●○○
SoundSmoothing™	✓	✓	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	✓	✓	✓
Auto EchoShield	✓	✓	Automatically reduces the effect of reverberations and sudden loud noises in reverberant environments.	✓		
Reverberant Room Program	✓	✓	Dedicated program for reverberant environments.	✓	✓	
Extended dynamic range	✓	✓	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓
Extended bandwidth	✓	✓	12 kHz available bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓		
Wearer Empowerment						
Signia Assistant	✓		Offers wearers instant 24/7 troubleshooting and customized sound adjustments via an intuitive AI-enabled chat. Provides detailed information in Connexx to address each client's fine-tuning needs more precisely.	✓	✓	✓
Signia App (iOS and Android)	✓	✓	Easy-to-use smartphone app enables remote control adjustments, audio streaming,** access to My WellBeing,** and support via Signia Assistant** and TeleCare.	✓	✓	✓
Adaptive Streaming Volume	✓		Speech intelligibility weighted adaptation of the streaming volume depending on the level of ambient noise.	✓	✓	✓
My WellBeing	✓		Accessed via the Signia App, wearers can get insights about how much they are using their hearing aids, their level of social interaction, their daily step count, and their overall activity level.	✓	✓	✓
Spatial Configurator	✓		Allows the wearer to adjust the direction of the microphone beam via the Signia App.	✓	✓	
TeleCare	✓	✓	Provides remote services including the TeleCare Dashboard, Remote Fine Tuning,** and CareChat (convenient text messages, phone calls, and video calls).	✓	✓	✓

*) Available for bilateral fittings with e2e wireless link
**) Requires Bluetooth connection, not available in Silk Charge&Go IX

Integrated Xperience feature overview

Features	Form Factor		Description	Performance levels		
	Pure	Silk		7 IX	5 IX	3 IX
Connectivity						
Direct streaming for iOS (MFi) & Android (ASHA)	✓		Hearing aids can be paired to Apple iOS devices (iPhones, iPads, iPods) and ASHA-supported Android devices for direct streaming of audio such as phone calls and music.	✓	✓	✓
HandsFree for iOS	✓		Enables truly hands-free calling. The wearer's voice is picked up by the hearing aid microphone and delivered in high quality to the conversation partner, while the call audio is streamed directly to the wearer's hearing aids.	✓	✓	✓
CallControl for iOS	✓		Calls can be answered, ended, and ignored with a simple press of the hearing aid button.	✓	✓	✓
Telecoil	✓		Available in Pure Charge&Go T IX only. Enables direction connection to audio loop systems and compatible telephone signals.	✓	✓	✓
StreamLine Mic compatibility	✓		Stream high quality stereo sound directly from any Bluetooth device (PCs, older Apple devices, non-ASHA Android devices, etc). Can also be used as a remote mic.	✓	✓	✓
StreamLine TV compatibility	✓		Stream high quality stereo sound directly from TV to hearing aids via Bluetooth.	✓	✓	✓
miniPocket compatibility	✓	✓	Optional keychain remote control.	✓	✓	✓
Tinnitus						
Notch Therapy	✓	✓	Signia's unique, patented tinnitus therapy which is proven especially effective against tonal tinnitus.	✓	✓	✓
Tinnitus Noise Therapy Manager	✓	✓	A choice of tinnitus treatments based on traditional noise therapy signals, amplified notch therapy, or a combination of the two approaches.	✓	✓	✓
Fitting Optimization						
Processing Channels / Gain & MPO handles	✓	✓	Number of processing channels/handles in Connexx	48/20	36/18	32/16
Frequency Compression	✓	✓	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓
CROS/BiCROS compatibility	✓	✓	A CROS device worn on an ear with unaidable hearing loss transmits sound to a hearing aid worn on the better ear.	✓	✓	✓
Hearing programs	✓	✓	Number of programs configurable in Connexx.	6	6	6
HD Music	✓	✓	Preset optimized music programs for enhanced enjoyment of non-streamed music.	3	3	1
Smart Optimizer and Data Logging	✓	✓	Tailored fine-tuning recommendations based on logging and analysis of the wearer's hearing aid use.	✓	✓	✓
Acclimatization Manager	✓	✓	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓
InSituGram	✓	✓	Allows quick and convenient measurement of hearing thresholds with the hearing aids in place.	✓	✓	✓
AutoFit™	✓	✓	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets, enabling significantly faster real-ear verification.	✓	✓	✓

*) Available for bilateral fittings with e2e wireless link

**) Requires Bluetooth connection, not available in Silk Charge&Go IX




Seamless hearing with discretion



Key features

-  Binaural OneMic Directionality
-  On-the-go rechargeability
-  Instant fit
-  Discreet look
-  Signia app
-  CROS compatibility



One push and you're done.
The EarWear 3.0 sleeves come in four different sizes and a closed and open build version suitable for most ear anatomies and hearing needs.



Silk IX Charger
With contactless charging, wearers simply place Silk Charge&Go IX in the charger without having to worry about aligning any charging contacts.

Silk Charge&Go IX

Performance levels

- 7IX
- 5IX
- 3IX
- sDemo

Comfortable and discreet, this tiny hearing aid provides the convenience of on-the-go rechargeability. It is packed with advanced sound features to support you in any environment including its unique Binaural OneMic Directionality, which offers unmatched conversation support in noisy environments for a hearing aid of its size.

Also available as a CROS Silk Charge&Go IX.

Click ITE	
Model	Click CIC
	
Performance	50/114 dB SPL
Battery	
Type	Lithium-ion power cell

Deliver discreet, advanced solutions



Key features

-  RealTime Conversation Enhancement
-  Augmented Focus
-  Own Voice Processing 2.0
-  Signia Assistant
-  Android & iPhone connectivity
-  Convenient rechargeability
-  Signia app
-  CROS compatibility
-  Available as a T-Coil option



Pure Charge&Go IX / Pure Charge&Go T IX

Performance levels

- 7IX
- 5IX
- 3IX
- sDemo

Make sure you never miss a word with RealTime Conversation Enhancement. Small, discreet, and conveniently rechargeable for all-day listening, with Bluetooth connectivity to smartphones. Pure Charge&Go T IX offers T-Coil and extended battery life.

Also available as a CROS Pure Charge&Go IX, or as Pure Charge&Go T IX with integrated T-Coil.

3 charging options
The Pure Portable Charger offers comfort and charging on the go. The Pure Dry&Clean Charger also offers UV cleaning and drying to keep the hearing aids in peak condition. And the small, reliable standard Pure Charger meets basic desktop charging needs.

RIC				
Model				
Receiver type	miniReceiver 3.0 S	miniReceiver 3.0 M	miniReceiver 3.0 P	miniReceiver 3.0 HP
   				
Performance	46/110 dB SPL	60/119 dB SPL	65/122 dB SPL	75/131 dB SPL
Battery				
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

Signia AX

Revolutionary split processing

Performance levels

- 7AX
- 5AX
- 3AX
- 2AX
- 1AX

Hearing isn't always easy. A big group of people talking at the same time, softly spoken talkers, too much background noise... Sometimes sounds blend together making it impossible to focus on what you want to hear.

Augmented Xperience changes the way your clients hear the world. Instead of them having to strain to discern speech from background noise, our new platform uses proprietary Augmented Focus™ technology to process the two separately and create a clear contrast. It then recombines them to deliver outstanding speech clarity in a fully immersive environment. For instance, 100% of participants in a recent study¹ reported excellent speech understanding in their home environment using this technology.



Augmented Xperience is available across all performance levels from 7 to 1 so you can offer all your customers outstanding speech clarity in an immersive environment

ITE



Insio Charge&Go AX

RIC



Pure Charge&Go AX



Pure Charge&Go T AX



Pure 312 AX

SLIM-RIC



Styletto AX

(1) Jensen, Pischel, Taylor, Schulte: Performance of Signia AX in At-Home Listening Situations (Signia White Paper, 2021)





Augmented Xperience feature overview

Our leading-edge Augmented Xperience platform is the first to offer outstanding speech clarity in an immersive environment, combined with connectivity for Android¹ and iOS devices plus artificial intelligence (AI) for live remote assistance.

Features	Information	Performance levels				
		7AX	5AX	3AX	2AX	1AX
Improved AX Soundscape Processing	Dynamic Soundscape Processing takes sound processing to a more precise level by including motion data for a more realistic sound adapted to everyday life.	●●●●●	●●●●●	●●●●●	●●●○○	●●●○○
Augmented Focus™	Separates the focus sound (speech) from surroundings, creating a clear contrast between the two, which makes it easier to follow conversations.	✓	✓	✓	✓	✓
Acoustic Sensors	The redesigned acoustic sensors can detect more variables of sound than ever before, providing more accurate information about the environment.	✓	✓	✓	✓	✓
e2e wireless 4.0	Control coupling, synchronization, binaural directionality, Spatial Speech Focus, CROS/ BiCROS.	✓	✓	✓	✓	✓
Motion Sensor	An integrated motion sensor detects whether the wearer is walking, running or sitting still. This motion data complements the acoustic data to optimize Dynamic Soundscape Processing 2.0.	✓	✓	✓		
OVP2.0 ³	OVP2.0 utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression. Available for RIC/BTE devices with e2e wireless in binaural fittings.	✓	✓	✓		
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearing situations.	●●●●●	●●●●○	●●●○○	●●●○○	●●○○○
Signal processing / Gain & MPO	Channels/handles in Connexx. For Signia AX, channels are per processing chain (Focus/Surroundings).	48/20	32/16	24/12	16/8	16/8
Hearing Programs	Number of programs configurable in Connexx.	6	6	6	4	4
Extended dynamic range	117 dB 24 bit processing for enhanced digital encoding of soft and loud sounds.	✓	✓	✓	✓	✓
Speech and noise management	Reduces listening effort using sophisticated sound processing.	✓	✓	✓	✓	✓
SoundSmoothing™	Reduces the annoyance of sudden loud impact sound such as when handling glassware.	✓	✓	✓	✓	
Feedback cancellation	High speed monitoring and control of feedback in individual processing channels.	✓	✓	✓	✓	✓
HD Music	Preset programs for enhanced enjoyment of non-streamed music.	3	3	1	1	
eWindScreen	Reduces the annoyance of wind noise when outdoors.	✓	✓	✓	✓	
Extended bandwidth	12 kHz bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓				
EchoShield	Dedicated program for reverberant environments.	✓	✓			
Speech Quality	Speech Quality delivers improved speech understanding and communication ease in noise.	●●●●●	●●●●○	●●●○○	●●●○○	●●●○○
Binaural Directionality	Directionality provides scientifically proven enhanced speech understanding in noise. Available for bilateral fittings with directional microphones and e2e wireless link.	✓	✓	✓		
Wireless CROS/BiCROS	Requires CROS transmitters.	✓	✓	✓	✓	✓
Frequency compression	Improves the detection of speech sounds for wearers with reduced high frequency hearing thresholds. Frequency compression is selectively activated based on audiological criteria.	✓	✓	✓	✓	✓

Features	Information	Performance levels				
		7AX	5AX	3AX	2AX	1AX
Spatial SpeechFocus	Depending on the direction of the dominant speech source, Spatial SpeechFocus steers the directional beam to the left, right, or behind the wearer automatically in a car situation, or in a dedicated "Strolling" program. Available for bilateral fittings with directional microphones and e2e wireless link.	✓ ¹⁾	✓ ²⁾			
Wearer Interaction	Technically inclined wearers may appreciate having extra control over their hearing aids. Signia Augmented Xperience hearing aids offer these wearers a comprehensive range of adjustment options via the Signia app.	●●●●●	●●●●●	●●●●●	●●●●●	●●●●●
Signia Assistant	Signia Assistant offers wearers instant 24/7 support and troubleshooting around the clock. It provides detailed information in Connexx to address each client's fine-tuning needs more precisely.	✓	✓	✓	✓	✓
Signia app (iOS and Android)	The new Signia app combines all functionalities of previous apps (touchControl, myControl, myHearing app). The functionality offered in the app is automatically customized depending on the connected hearing aid.	✓	✓	✓	✓	✓
Adaptive Streaming Volume	Speech intelligibility weighted adaptation of the streaming volume depending on ambient noise.	✓	✓	✓	✓	✓
Spatial Configurator	Spatial Configurator is part of the Signia app and allows the wearer to adjust the direction of the microphone beam.	✓	✓			
Direct streaming for iOS	Signia Augmented Xperience hearing aids can be connected directly to Apple devices for phone calls and audio streaming – without additional accessories.	✓	✓	✓	✓	✓
Direct streaming for Android	Direct audio streaming from ASHA compatible Android devices to Bluetooth enabled Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
StreamLine TV accessory	Stream high quality stereo sound directly from a TV to Bluetooth equipped Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
Call Control/ Tap Control + HandsFree for iOS	Call Control/Tap Control: accept and end call using a hearing aid Handsfree for iOS: hands-free mobile phone calls – even in noisy environments, combine phone calls with other activities or to listen to music or podcasts via your hearing aids	✓	✓	✓		
StreamLine Mic accessory	Stream high quality stereo sound directly from any Bluetooth device (Android devices, older Apple devices, other smartphones, PCs, laptops, etc.) to Bluetooth equipped Augmented Xperience hearing aids.	✓	✓	✓	✓	✓
Tinnitus	A choice of tinnitus treatments based on amplified notch therapy or traditional noise therapy signals. It is even possible to combine the two approaches.	✓	✓	✓	✓	
Notched Amplification Therapy	Innovative tinnitus therapy delivered via everyday amplification.	✓	✓	✓	✓	
Smart Optimizer and Data Logging	Smart Optimizer (found under Data Logging) provides tailored fine-tuning recommendations based on an analysis of the wearer's usage of the hearing aid.	✓	✓	✓	✓	✓
Acclimatization manager	Allows the wearer to benefit from enhanced audibility in the longer term, while enjoying optimal device acceptance when initially fitted.	✓	✓	✓	✓	✓
InSituGram	Allows quick and convenient screening audiograms and checking of thresholds while the hearing aids are worn.	✓	✓	✓	✓	✓
AutoFit	Automatically provides an optimized real-ear insertion gain (REIG) setting without the need to manually match targets. The new IMC2 standard expands the range of REM equipment which supports AutoFit with Connexx and Noah 4.5.1 or higher.	✓	✓	✓	✓	✓
TeleCare	TeleCare provides remote services including additional tools to follow-up with patients. In addition, user engagement and autonomy can be increased via the Signia app.	✓	✓	✓	✓	✓
Remote Services	Includes TeleCare Dashboard, Remote Fine Tuning, CareChat (incl. text messages, phone calls and video calls), virtual follow-up.	✓	✓	✓	✓	✓

1) In Universal program: Automatic mic beam to the right / left / rear directions.
2) In Universal program: Automatic mic beam to the rear, surround mode for side speech. In manual programs, microphone beam to the right/left/rear directions.
3) Not available for Insio Charge&Go AX



Designed for brilliant hearing and looks



Key features

- Augmented Focus
- Own Voice Processing 2.0
- Signia Assistant
- Android & iPhone connectivity
- Recharge on-the-go
- CROS compatibility
- Award-winning design*

Performance levels

- 7AX
- 5AX
- 3AX
- 2AX
- 1AX
- sDemo

Styletto AX

Our new Styletto AX rechargeable hearing aid looks as amazing as its sounds. Its award-winning SLIM-RIC design* is smaller than before for an even more discreet fit. Built on our leading-edge Augmented Xperience platform, Styletto AX uses our proprietary Augmented Focus™ technology to split speech from surrounding sounds, process them separately, and recombine them to deliver outstanding speech clarity in an immersive environment.

Also available as a CROS Styletto AX.



Portable charging options
Styletto AX offers the choice of a portable charger with Qi technology – wearers simply place the charger on a charging pad instead of having to plug it into the mains. Alternatively, they can opt for a Dry&Clean charger, which simultaneously charges the hearing aids and dries and cleans them via UV light.

SLIM-RIC

Model			
Receiver type	miniReceiver 3.0 S	miniReceiver 3.0 M	miniReceiver 3.0 P
Performance			
Matrix (2 ccm)	46/110 dB SPL	60/119 dB SPL	65/122 dB SPL
Battery			
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

* Red Dot Design Award 2019 for Styletto Connect, iF Design Award 2019 for Styletto



Immersive sound in our smallest rechargeable RIC



Key features

- Augmented Focus
- Own Voice Processing 2.0
- Signia Assistant
- Android & iPhone connectivity
- Convenient rechargeability
- CROS compatibility
- Available as T-Coil Option

Performance levels

- 7AX
- 5AX
- 3AX
- 2AX
- 1AX
- sDemo

Pure Charge&Go AX / Pure Charge&Go T AX

Now available in all performance levels, Pure Charge&Go AX is ideal for experienced and first-time hearing aid wearers who are looking for a sleek aesthetic device that is fully featured. It combines an attractive design with proven rechargeability, state-of-the-art Bluetooth connectivity for Android and iPhone, and the landmark AI of the Signia Assistant.

Also available as a CROS Pure Charge&Go AX, or as Pure Charge&Go T AX with integrated T-Coil.



3 charging options
The Pure Portable Charger offers comfort and charging on the go. The Pure Dry&Clean Charger also offers UV cleaning and drying to keep the hearing aids in peak condition. And the small, reliable standard Pure Charger meets basic desktop charging needs.

RIC

Model				
Receiver type	miniReceiver 3.0 S	miniReceiver 3.0 M	miniReceiver 3.0 P	miniReceiver 3.0 HP
Performance				
Matrix (2 ccm)	46/110 dB SPL	60/119 dB SPL	65/122 dB SPL	75/131 dB SPL
Battery				
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell



The first custom hearing aids with contactless charging and Bluetooth



Performance levels

- 7AX
- 5AX
- 3AX
- 2AX
- 1AX

Insio Charge&Go AX

Our new rechargeable Insio Charge&Go AX ITC and ITE models redefine what custom hearing aids can offer experienced wearers. The smaller, more discreet design offers a higher fitting rate than previous devices and all-day comfort, while our leading-edge Augmented Xperience platform delivers outstanding speech clarity in an immersive environment. The custom hearing aids even include a motion sensor to adapt their sound optimally to the wearer’s situation by detecting if they are moving or stationary.

- Key features
- Augmented Focus
 - Signia Assistant
 - Android & iPhone connectivity
 - Contactless charging
 - Custom fit
 - Award-winning design*



Insio Charger
These are the first custom hearing aids with contactless charging. Wearers simply place Insio Charge&Go AX in the Insio Charger without having to worry about any exact charging contact.

	ITC/HS ITE	ITC/HS ITE	ITC/HS
Model			
Performance	124/65 dB SPL	118/55 dB SPL	113/50 dB SPL
Battery	Lithium-ion power cell		
Type			

* CES Innovation awards 2022 honoree and iF Design Award 2022 for Insio Charge&Go AX



The RIC for outstanding speech clarity



Performance levels

- 7AX
- 5AX
- 3AX
- 2AX
- 1AX
- sDemo

Pure 312 AX

Equipped with all the benefits of our Augmented Xperience platform for outstanding speech clarity, Pure 312 AX and CROS Pure 312 AX attract wearers with their fresh, elegantly slim design. Despite their discreet form, these tiny hearing aids provide fully-featured Bluetooth connectivity for maximum enjoyment of music, calls and TV from Android and iOS devices, plus all the benefits of our Signia Assistant. Pure 312 AX’s Zn-Air battery delivers a runtime of around 4 days.*

- Key features
- Augmented Focus
 - Own Voice Processing 2.0
 - Signia Assistant
 - Android & iPhone connectivity
 - CROS compatibility


	RIC			
Model	miniReceiver 3.0 S	miniReceiver 3.0 M	miniReceiver 3.0 P	miniReceiver 3.0 HP
Receiver type				
Performance	44/108 dB SPL	58/116 dB SPL	62/119 dB SPL	74/128 dB SPL
Matrix (2 ccm)				
Battery	312	312	312	312
Size	1.5mA	1.7mA	1.7mA	1.7mA
Power consumption				

*Based on 17 hours wearing per day including 5 hours of streaming

Signia Xperience

Empowering people to hear what matters

Performance levels



How do you convince people with hearing loss who have not yet decided to wear hearing aids to do so? At Signia, we have dedicated ourselves to this challenge by finding ways for our solutions to help their wearers Be Brilliant. People with hearing loss do not want to feel old because they aren’t old. They want hearing solutions that match their active lifestyle and can-do spirit. Our Signia X hearing aids offer leading audiological performance in iconic form factors. They combine our proven Signia Xperience YourSound technology built on revolutionary acoustic-motion sensors with Li-ion rechargeability and Bluetooth connectivity so your clients always hear what matters, with ease.



Xperience

Signia Xperience is available across all performance levels from 7 to 1, so you can offer all your customers outstanding hearing performance.

Click ITE

RIC



Active / Active Pro



Silk X



Pure Charge&Go X



Pure 312 X

SLIM-RIC

BTE



Styletto X



Motion Charge&Go X



Motion Charge&Go P X



Motion Charge&Go SP X

Signia Xperience feature overview

Our landmark platform Signia Xperience gives hearing aids a new sense to truly understand what the wearer needs to hear even in noise, all combined with advanced connectivity and leading remote assistance.

Features	Information	Performance levels						
		7X	5X	3X	2X	1X	Active Pro	Active
YourSound Technology	Redesigned acoustic sensors plus an integrated motion sensor provide detailed information about where the wearers are and what they are doing. Dynamic Soundscape Processing and Own Voice Processing (OVP™) process this information delivering the most clear and personalized sound.	●●●●●	●●●●●	●●●●●	●●●○○	●●●○○	●●●●●	●●●○○
Acoustic Sensors	Detects more variables of sound than ever before, providing more accurate information about the environment.	✓	✓	✓	✓	✓	✓	✓
Dynamic Soundscape Processing	Steers Sound Clarity and Sound Quality features for natural sound and speech in every situation from any direction, even when moving.	●●●●●	●●●●●	●●●●●	●●●○○	●●●○○	●●●○○	●●●○○
Motion Sensor ^A	Detects whether the wearer is walking, running or not. This information is used for Dynamic Soundscape Processing.	✓	✓	✓				
OVP™ (Own Voice Processing) ^A	Utilizes real-time recognition of the wearer’s voice to deliver a natural own-voice impression.	●●●●●	●●●●●	●●●●●				
Signal Processing / Gain & MPO	Channels/handles in Connexx.	48/20 ^B	32/16	24/12	16/8	16/8	48/20	16/8
Signia Assistant ^A	Offers the wearers instant 24/7 support and troubleshooting around the clock. It provides detailed information in Connexx to address each client’s fine-tuning needs more precisely.	✓	✓	✓	✓	✓	✓	✓
TeleCare	Provides virtual appointments, remote fine-tuning, progress updates, hearing lessons and daily satisfaction ratings.	✓	✓	✓	✓	✓	✓	✓
Tinnitus	A choice of tinnitus treatments based on amplified notch therapy or traditional noise therapy signals. It is even possible to combine the two approaches.	✓	✓	✓	✓		✓	
eWindScreen™	Reduces the annoyance of wind noise when outdoors, improving sound comfort. eWindScreen Binaural requires an e2e wireless link.	✓ Binaural ^A Auto in Universal program	✓ Binaural ^A manual	✓ Monaural	✓ Monaural			✓ Monaural
Binaural OneMic Directionality 2.0	Delivers a highly sophisticated speech-in-noise performance. Available for Silk X.	✓	✓	✓				
Spatial SpeechFocus	Detects the loudest voice in the environment and adjusts the directional beam towards that voice (i.e. front, side, or backward facing). Available for bilateral fittings with directional microphones and e2e wireless link.	✓ front / back / right / left ^C	✓ front / back / omni ^D				✓ front / back / right / left ^C	
Extended Bandwidth	12 kHz bandwidth for enhanced processing of high-frequency speech and environmental sounds.	✓						✓
EchoShield ^E	Dedicated program for reverberant environments.	✓						✓

●●●●● = highest feature performance
●●●● = mid-range feature performance
●●● = base feature performance

A) Not available with Silk X
B) For Motion Charge&Go SP X: 34/16
C) In Universal program: automatic mic beam to the right / left / rear directions
D) In Universal program: automatic mic beam to the rear, surround mode for side speech. In manual programs, microphone beam to the right / left / rear directions

E) Study conducted at the University of Northern Colorado, 2015, examined the effectiveness of new Signia features by collecting and analyzing ongoing EEG data while subjects performed speech testing. For both Signia features SpeechMaster and EchoShield, the objective brain behavior measures revealed a significant reduction in listening effort when the feature was activated.




Prescription audio that performs

Key features

- Iconic earbud design bridging the gap to high tech consumer products
- Signia Xperience technology for brilliant sound and speech understanding
- Remote control
- TeleCare

- Signia app ready
- Instant fit
- Equipped with Bluetooth
- Recharge on-the-go
- Dynamic Soundscape Processing



Portable Charger
Signia Active Pro comes in a pocket-sized Qi wireless charging case that lets wearers put the buds in and out of their ears when they need them.

Active / Active Pro

Performance levels

- Active Pro
- Active

Our game-changing Signia Active Pro takes hearing aids to that next level for you: Prescription audio that enhances human performance and enables wearers to Be Brilliant. Signia Active can be put in and taken out of the ears as easily as any earbud. Modern-looking, easy to use and fast to fit, these handy devices mark a whole new category of hearing solutions – for everyone who doesn’t want to wear hearing aids but wants to enjoy hearing everything that’s important to them.

Click ITE	
Model	
Performance	
Matrix (2 ccm)	42/118 dB SPL
Battery	
Type	Lithium-ion power cell

Experience the full potential of BTEs

Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World’s first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV

- Signia app ready
- Smaller Design
- Equipped with Bluetooth
- Rechargeable
- Own Voice Processing



Inductive Charger II
This charger comes with a charger lid to protect the hearing aids, a fast charging option, and a drying function. Its smart design accommodates custom molds.

Motion Charge&Go X

Comparable in its small size to a 312- device, Motion Charge&Go X is the ultimate BTE combining sleek design with rechargeability, Bluetooth, and the Signia Assistant. It is flexible to fit, easy to use and comfortable to wear, making this an easy choice for those looking for convenience and robustness – all in one.

Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

BTE			
Model	Earhook	ThinTube 3.0	ThinTube 3.0 P
Performance			
Matrix (2 ccm)	60/133 dB SPL	60/125 dB SPL	63/126 dB SPL
Battery			
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

Hear what really matters



Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Motion Charge&Go P X

Smaller than its predecessor without compromising on power, Motion Charge&Go P X combines Signia Xperience hearing technology with the complete convenience of Li-ion recharging, making annoying battery changes a thing of the past. This powerful hearing solution for every situation with T-Coil ensures that whatever you spend your day doing, you will be able to perform at your best.

Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV



Signia app ready



Smaller Design



Equipped with Bluetooth



Rechargeable



Integrated T-Coil



Own Voice Processing



Charging Options
The consumer can select between two charger options, a standard charger or a Dry&Clean charger.

BTE			
Model	Earhook	ThinTube 3.0	ThinTube 3.0 P
			
Performance			
Matrix (2 ccm)	77/135 dB SPL	66/130 dB SPL	70/131 dB SPL
Battery			
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

Superpower your hearing



Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Motion Charge&Go SP X

This is the world's first rechargeable Super Power BTE. Closing the gap between Super Power and super features, the innovative Motion Charge&Go SP X ultra-power device is the complete package for those who want to Be Brilliant despite severe to profound hearing loss. With Bluetooth connectivity, Li-ion rechargeability and T-Coil, it offers a world of new hearing and handling possibilities not familiar to those with conventional SP devices.

Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV



Signia app ready



Smaller Design



Equipped with Bluetooth



Rechargeable



Integrated T-Coil



Own Voice Processing



Charging Options
The consumer can select between two charger options, a standard charger or a Dry&Clean charger.

BTE	
Model	Motion Charge&Go SP X
	
Performance	
Matrix (2 ccm)	82/140 dB SPL
Battery	
Type	Lithium-ion power cell

Hearing without compromise

- Performance levels
- 7X

5X

3X

2X

1X

sDemo

Styletto X

Styletto X is the complete package hearing aid ideal for anyone who doesn’t want to compromise on their active lifestyle just because of hearing loss. The world’s first SLIM-RIC combines brilliant, award-winning design* and portable rechargeability with exchangeable receivers, the clear sound and speech understanding of our revolutionary Signia Xperience platform, and Bluetooth connectivity. This unique form factor can significantly increase your conversion of people with hearing loss into satisfied clients. When consumers were given a selection of hearing aids to choose from, Styletto increased conversion by 18%¹.

SLIM-RIC			
Model			
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P
			
Performance			
Matrix (2 ccm)	45/108 dB SPL	60/119 dB SPL	70 /124 dB SPL
Battery			
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell

(1) Source: Hakvoort, Burton: Increasing Style, Reducing Stigma: The Styletto Solution (Signia Backgrounder, 2018)

* Red Dot Design Award 2019 for Styletto Connect, iF Design Award 2019 for Styletto

- Key features
-  Signia app ready

 Recharge on-the-go

 Equipped with Bluetooth

 Own Voice Processing



Portable Charger
The world’s first hearing aid charger with Qi technology – wearers simply place the charger on a charging pad instead of having to plug it into the mains.

Discretion guaranteed

- Performance levels
- 7X

5X

3X

2X

1X

sDemo

Silk X

Silk X is the ideal solution for those who value discretion above all. The only ready-to-wear CIC available on the Signia Xperience platform delivers clear speech understanding in an ultra-discreet size with instant fit. It is tailor-made without the wait: Your clients can walk in, experience Silk X and walk out the door with a new CIC without having to wait for a custom-built device. Thanks to its flexible silicone Click Sleeves, Silk X sits securely in the ear for a brilliant fit on the spot.


Click ITE	
Model	Click CIC
	
Performance	
Matrix (2 ccm)	50/114 dB SPL
Battery	
Size	10
Power consumption	1.2 mA

- Key features
-  Signia app ready

 Compatible with CROS device for people with unaidable hearing loss on one side

 Instant fit

 Dynamic Soundscape Processing



One click and you’re done.
The sleeves come in four different sizes and a closed and open build version suitable for most ear anatomies and hearing needs.







Supreme style and streaming with charging-on-the-go



Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming with StreamLine Mic and StreamLine TV



Inductive Charger II
This charger comes with a charger lid to protect the hearing aids, a fast charging option, and a drying function. Its smart design accommodates custom molds.

Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Pure Charge&Go X

Pure Charge&Go X combines all the advantages of Signia Xperience with Li-ion rechargeability and fully-featured Bluetooth connectivity. Its 20% extra battery capacity provides full flexibility and convenience for higher wearer satisfaction, with streaming of music, calls and TV in high-quality sound. All in a remarkably small design. Simply Charge&Go!

RIC

Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Performance				
Matrix (2 ccm)	45/108 dB SPL	60/119 dB SPL	70/124 dB SPL	75/130 dB SPL
Battery				
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell







Fully versatile, fully connected



Key features

- Most natural own voice and best speech understanding in noise, even when the wearer is in motion
- World's first integrated acoustic-motion sensors
- Long-lasting high quality audio streaming
- TeleCare
- Remote Control
- Connectivity & Streaming



Performance levels

- 7X
- 5X
- 3X
- 2X
- 1X
- sDemo

Pure 312 X

Equipped with all the benefits of the Signia Xperience platform, Pure 312 X attracts wearers with its elegantly slim design. Despite its discreet form, Pure 312 X provides fully-featured Bluetooth connectivity for maximum enjoyment of music, calls and TV. An integrated T-Coil is also available.

RIC

Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Performance				
Matrix (2 ccm)	45/108 dB SPL	60/119dB SPL	70/124 dB SPL	75/130 dB SPL
Battery				
Size	312	312	312	312
Power consumption	1.7 mA	1.9 mA	1.8 mA	1.8 mA

Signia Nx

Outstanding
soundscape
processing

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx

The sophisticated chip technology of Signia Nx packs even more performance into even less space. Its leading technology enables advanced connectivity in the smallest sized hearing aid with Ultra HD e2e and Bluetooth.



Signia Nx

Motion 13 Nx, 13P Nx and Insio Nx cover all performance ranges from performance levels 7 to 1. That means more choice for your customers, more opportunities for you.

Signia Nx is available across all performance levels from 7 to 1, so you can offer all your customers outstanding hearing performance.

BTE



Motion 13 Nx



Motion 13P Nx

ITE



Insio Nx

Signia Nx feature overview

Our game-changing platform Signia Nx delivers the most natural own voice and speech understanding in noise, while combining the most advanced connectivity systems.

Features	Information	Performance levels				
Audiology		7Nx	5Nx	3Nx	2Nx	1Nx
OVP™ (Own Voice Processing)	Own Voice Processing (OVP™) utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression. Available for SLIM-RIC/RIC/BTE devices with Ultra HD e2e.	●●●●	●●●●	●●●●	-	-
3D Classifier	Unlike conventional classification systems which rely on acoustic data, 3D Classifier analyzes multi-dimensional information from the wearer's environment and communication behavior. Information on the wearer's voice activity, motion state and acoustic data are used to tailor the hearing aids' response to the wearer's needs. Motion detection available for Bluetooth-enabled devices.	●●●●	●●●●	●●●	-	-
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearing situations. Exceptional clarity is achieved with features such as: HD Spatial, extended dynamic range, extended bandwidth, EchoShield, HD Music, eWindScreen™, speech and noise management and feedback cancellation.	●●●●	●●●●	●●●	●●	●●
Speech Quality	Speech Quality delivers improved understanding of other people's voices in noise by using an exclusive technology called Ultra HD e2e. Features which utilize Ultra HD e2e include narrow directionality, Spatial SpeechFocus™ and TwinPhone.	●●●●	●●●●	●●●	●●	●●
Connectivity						
Direct Streaming	Signia Nx hearing aids can be connected directly to Apple devices for phone calls and audio streaming – without additional accessories. With the StreamLine Mic accessory, direct streaming is possible also from Android devices and other Bluetooth mobile devices. With the StreamLine TV accessory, high-quality stereo sound can be streamed directly to Nx hearing aids. Available for Bluetooth-equipped devices.	✓	✓	✓	✓	✓
You Adjust	Technically inclined wearers often appreciate having extra control over their hearing aids. Signia Nx hearing aids offer these wearers a comprehensive range of adjustment options via Signia app.	●●●●	●●●●	●●●●	●●●●	●●●●
Remote hearing care						
TeleCare	TeleCare provides remote services including additional tools to follow-up the patients within the whole journey, beginning with pre-appointment and even continuing after purchase. It contributes to a more successful trial phase and ensures higher flexibility. In addition, the user engagement and autonomy during this whole journey will be increased via the Signia app.	✓	✓	✓	✓	✓

●●●●● = highest feature performance



Direct streaming with the most versatile fitting options



Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx
- sDemo

Motion 13 Nx

Available in a small model and a sleek P model, the incredibly versatile Motion 13 Nx BTE is ideal for wearers with individual fitting requirements and moderate to severe or even profound hearing loss. It harmonizes the wearer’s voice and surrounding sounds, providing the most natural own voice and best speech understanding in noise. Motion 13 Nx and Motion 13P Nx provide full Bluetooth connectivity and Signia TeleCare with an optional modular T-Coil integrated in the battery door.

Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- For individual handling and fitting requirements
- TeleCare
- Remote Control
- Connectivity & Streaming

Signia app ready

Equipped with Bluetooth

Standard battery door can be exchanged for one with an integrated T-Coil

Compatible with CROS device for people with unaidable hearing loss on one side

Own Voice Processing

Ultra HD E2E



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Motion 13 Nx hearing aids.

StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

BTE		
Model	Motion 13 Nx	Motion 13P Nx
Performance	ThinTube 56/126 dB SPL Earhook 63/134 dB SPL	ThinTube 65/130 dB SPL Earhook 77/135 dB SPL
Battery		
Size	13	13
Power consumption	1.4 mA	1.7 – 2.0 mA



The most discreet natural hearing experience with full connectivity

Key features

- Best speech understanding in noise and a most natural sound thanks to Nx chip
- Small size up to high power matrix
- Full wireless and remote control options for all models
- TeleCare
- Remote Control
- Connectivity & Streaming*

Signia app ready

Equipped with Bluetooth

Compatible with CROS device for people with unaidable hearing loss on one side

Ultra HD E2E




StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Bluetooth-enabled Insio Nx ITC and ITE hearing aids.

StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any Bluetooth device for Bluetooth-enabled Insio Nx ITC/ITEs.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- 2Nx
- 1Nx



Insio Nx

Signia’s Insio™ Nx ITC and ITE models provide top Bluetooth connectivity with best-in-class sound quality. Tailor-made for the most discerning wearers, Insio Nx delivers a more natural sound with maximum discretion. Equipped with Ultra HD e2e for superior wireless binaural processing, including as a small CIC with a single microphone, Insio offers clear sound and speech understanding even in noisy environments.

ITE				
Model	IIC	CIC	ITC	ITE
Performance				
Matrix (2 ccm)	50/113 dB SPL ¹	S 50/113 dB SPL ¹ M 55/118 dB SPL ¹ P 65/124 dB SPL ¹	S 50/113 dB SPL ¹ M 55/118 dB SPL ¹ P 65/124 dB SPL ¹	M 55/118 dB SPL ¹ P 65/124 dB SPL ¹
Battery				
Size	10	10	312	312
Power consumption	1.3 mA ¹	1.3 mA ¹	1.3 mA ¹	1.3 mA ¹

(1) Max. values


* Available with Bluetooth-enabled Insio Nx ITC and ITE hearing aids

Essential

Proven technology,
outstanding
comfort

Whether it’s a receiver-in-canal, RIC or behind-the-ear, BTE or one of the in-the-ear, ITE models, the Intuis 3 family has a solution for you. Suitable for mild to severe hearing loss, all hearing aids provide you with a reliable, clear and comfortable listening experience and a maximum of discretion with remote control options via the Signia app or the miniPocket™.

Hearing what matters has never been more important than in today’s busy world, where a missed word can mean a missed opportunity. Signia Intuis 4 hearing aids are packed with the power, modern functionality and support you need to truly shine. Suitable for mild to moderate hearing loss, Intuis M 4’s state-of-the-art sensors and microphones help you to enjoy brilliant sound and outstandingly clear speech, so you can easily follow conversations even in challenging situations with a lot of background noise.



Essential segment feature overview

Proven technology, outstanding comfort: In a broad range of hearing aid models, the Essential line offers a reliable hearing experience, a high level of wearing comfort and easy handling for almost every customer.

Features	Intuis 4			Intuis 3
	Essential			
	4.2	4.1	4.0	
Augmented Focus	●	●	●	-
Acoustic sensor	●	●	●	-
Motion sensor	-	-	-	-
Frequency channels / compression controls	16/8	16/8	12/6	12/6
Number of programs	4	4	4	4
Extended dynamic range	●	●	●	-
Speech and noise management	●	●	●	●
SoundSmoothing	●	-	-	-
Feedback cancellation	●	●	●	●
HD Music (presets)	1	-	-	-
eWindScreen	●	-	-	-
Directional microphones	Adaptive	Adaptive	automatic	automatic ^B
Wireless CROS/BiCROS	●	●	-	-
Frequency compression	●	●	-	-
e2e wireless	●	●	●	●
Signia Assistant	●	●	-	-
Signia app	●	●	●	●
Adaptive streaming volume	●	●	-	-
Direct streaming Android devices (ASHA)	●	●	-	-
Direct streaming iOS	●	●	-	-
HandsFree for iOS	-	-	-	-
Tinnitus therapy	●	-	-	-
Smart optimizer and Data logging	●	●	●	●
Acclimatization manager	●	●	-	-
InSituGram	●	●	●	●
AutoFit	●	●	●	●
miniPocket	●	●	●	●
StreamLine TV and Mic	●	●	-	-

B) Only BTE

Small and powerful



- Key features
- For mild to severe hearing loss
 - Small design for a high level of wearing comfort
 - TeleCare
 - Remote Control
 - 6 channel MPO for the most advanced fitting options and more individual adjustments
- 

Signia app ready
- 

T-Coil



Signia app and miniPocket
To enhance its discretion, Intuis 3 RIC312 is remotely controllable via the Signia app and the miniPocket.

Intuis 3 RIC312

Performance levels

Essential

The remotely controllable Intuis 3 RIC312 is an attractive and convenient solution for mild to severe hearing loss. Thanks to its automatic directionality and a variety of noise reduction features, it offers a reliable hearing experience and improved speech intelligibility even in noisy environments. TeleCare lets your customers enjoy live remote support and fast, smooth trials.

RIC				
Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
				
Performance				
Matrix (2 ccm)	45/108 dB SPL	60/119 dB SPL	70/124 dB SPL	75/130 dB SPL
Battery				
Size	312	312	312	312
Power consumption	0.9 mA	1.0 mA	1.0 mA	1.1 mA

The latest technology for more of your customers



- Key features
- For all levels of hearing loss
 - Brilliant sound
 - Bluetooth streaming¹
 - Easy-to-use controls
 - TeleCare
- 

Signia app ready
- 

Equipped with Bluetooth
- 

Integrated T-Coil
- 

Ultra HD E2E
- 

Dynamic Soundscape Processing



Choose from a range of colors

Intuis 4 BTEs

Performance levels

PL 4.2

PL 4.1

PL 4.0

Intuis 4 BTE hearing aids feature the latest audiological technology, providing brilliant sound and speech clarity, while enhanced connectivity and Bluetooth streaming provide a hearing experience comparable to that of our most advanced devices.¹ Suitable for all levels of hearing loss, the ergonomic design and easy-to-use controls of Intuis 4 make daily operation straightforward and comfortable. The Signia app provides a suite of controls and communication options - so even more of your customers can take full advantage of the latest technology.

BTE			
Model			
	Intuis M 4	Intuis P 4	Intuis SP 4
			
Performance			
Matrix (2 ccm)	ThinTube 3.0: 61/125 dB SPL Earhook: 60/133 dB SPL	ThinTube 3.0: 66/129 dB SPL Earhook: 77/135 dB SPL	Earhook: 82/140 dB SPL
Battery			
Size	13	13	675

¹⁾ Bluetooth connectivity is available at PL 4.2 and PL 4.1

Proven technology, outstanding comfort

- Key features
- For all levels of hearing loss
- Ergonomic design for a high level of wearing comfort
- Remotely controllable for a maximum of convenience and discretion
- TeleCare
- Remote Control
- 6 channel MPO for the most advanced fitting options and more individual adjustments

Signia app ready

T-Coil (except Intuis 3 S)



Intuis 3 BTEs

Performance levels

Essential

All Signia’s Intuis 3 BTE models offer proven technology with 6 channel MPO for the most advanced fitting options and for more individual adjustments. Suitable for all levels of hearing loss, the easy to handle Intuis 3 BTEs are equipped with automatic directional microphones for a precise hearing experience, offer remote control options, and are TeleCare ready for smoother, faster trials and higher customer satisfaction.

BTE				
Model	Intuis 3 S	Intuis 3 M	Intuis 3 P	Intuis 3 SP
				
Performance	LifeTube 45/124 dB SPL Earhook 55/124 dB SPL	LifeTube 53/125 dB SPL Earhook 60/130 dB SPL	Earhook 70/134 dB SPL	Earhook 80/140 dB SPL
Battery				
Size	312	13	13	675
Power consumption	0.9 mA	1.0 mA	1.4 mA	2.4 mA

The perfect match

- Key features
- For mild to severe hearing loss
- 6 channel MPO for the most advanced fitting options and more individual adjustments
- Remotely controllable for a maximum of discretion and convenience
- TeleCare
- Remote Control

Signia app ready

Optional T-Coil



Intuis 3 Click ITEs and ITEs

Performance levels

Essential

With Intuis 3 Click ITEs and tailor-made ITEs, you can offer your customers the full range of in-the-ear models: No matter if they wish for a solution that fits right away or prefer a customized model, Intuis 3 is the perfect match. Equipped with proven technology, all models offer a precise hearing experience, are remotely controllable for a maximum of discretion and TeleCare ready for successful trials.

Click ITE/ITE		
Model	Intuis 3 Click CIC, Click ITC	Intuis 3 CIC, ITC, ITE
		
Performance	Click CIC 50/113 dB SPL Click ITC 55/116 dB SPL	CIC 65/124 dB SPL ¹ ITC 65/124 dB SPL ¹ ITE 65/124 dB SPL ¹
Battery		
Size	Click CIC 10, Click ITC 312	CIC 10, ITC 312, ITE 13
Power consumption	1.1 mA	1.1 mA ¹

(1) Max. values

Basic

Help the world hear better

The Basic segment of hearing aids offers fitting opportunities that enable you to bring the promise of better hearing to people who would otherwise have no hope of a solution. Easy to learn and easy to fit, even in mobile situations, you can take better hearing to where it’s needed – even to outlying areas. This range of BTE and Click ITE devices is suitable for all levels of hearing loss, sizes of ear canal and individual preference.

Basic

Basic segment feature overview


Available in a wide range of fittings and suitable for all levels of hearing loss from mild to profound, the Basic range offers a good hearing experience for every customer.

Features	Prompt						
	S	P	SP	Click ITEs	CIC	ITC	ITE
Matrix	LifeTube 45/124 dB Earhook 55/124 dB	70/134 dB	80/140 dB	Click CIC 50/113 dB Click ITC 55/116 dB	65/124 dB ¹	65/124 dB ¹	65/124 dB ¹
Channels/Controls	8/4	8/4	8/4	8/4	8/4	8/4	8/4
Programs	4	4	4	Click CIC: 4 ² Click ITC: 4	4	4	4
Directional microphones	●	●	●	-	-	-	-
T-Coil	-	●	●	-	-	○	○
Audio input	-	○	○	-	-	-	-
Noise reduction	situation based	situation based	situation based	situation based	situation based	situation based	situation based
SoundBalance™	-	-	-	-	-	-	-
Feedback cancellation	●	●	●	●	●	●	●

BTE



Click ITE



ITE



Features	Run				Fun		Fast
	S	P	SP	Click ITEs	P	SP	P
Matrix	LifeTube 45/124 dB SPL Earhook 55/124 dB SPL	70/132 dB SPL	79/140 dB SPL	Click CIC 50/113 dB SPL Click ITC 55/116 dB SPL	70/132 dB SPL	79/140 dB SPL	62/130 dB SPL
Channels	8	8	8	8	6	6	4
Programs	6 preconfigured sound profiles	6 preconfigured sound profiles	6 preconfigured sound profiles	6 preconfigured sound profiles	4 preconfigured sound profiles	4 preconfigured sound profiles	3 preconfigured sound profiles
Directional microphones	●	●	-	-	-	-	-
T-Coil	-	-	-	-	-	-	-
Audio input	-	-	-	-	-	-	-
Noise reduction	situation based	situation based	situation based	situation based	always on	always on	always on
SoundBalance™	-	-	-	-	-	-	-
Feedback cancellation	●	●	●	●	●	●	●


(1) Max. values
(2) 1 program + 4 accessible via app



Make your customers the center of attention

Key features

- A full range of BTE, Click ITE and ITE models
- For all levels of hearing loss
- Remotely controllable for a maximum of convenience and discretion
- TeleCare
- Remote Control



Basic

The Prompt family

Signia’s Prompt family offers outstanding wearing comfort and ease of use for every form of hearing loss. All models offer remote control options via the Signia app and the miniPocket™ and are a convenient solution for almost all your customers.

BTE			
Model	Prompt S	Prompt P	Prompt SP
Performance			
Matrix (2 ccm)	ThinTube 45/124 dB SPL Earhook 55/124 dB SPL	Earhook 70/134 dB SPL	Earhook 80/140 dB SPL
Battery			
Size	312	13	675
Power consumption	0.9 mA	1.4 mA	2.4 mA

Click ITE/ITE

Model		Prompt Click CIC, Click ITC	Prompt CIC, ITC, ITE
Performance			
Matrix (2 ccm)	Click CIC 50/113 dB SPL Click ITC 55/116 dB SPL		CIC 65/124 dB SPL ¹ ITC 65/124 dB SPL ¹ ITE 65/124 dB SPL ¹
Battery			
Size	Click CIC 10, Click ITC 312		CIC 10, ITC 312/10, ITE 13/312
Power consumption	1.1 mA		1.1 mA ¹


(1) Max. values



Take hearing to a new level

Key features

- Quick and easy fittings everywhere via the MobileFitting app
- For mild to profound hearing loss
- Ergonomic shape for a high level of wearing comfort
- TeleCare
- Remote Control



Basic

Run

Easy to fit in just a few minutes with the convenient MobileFitting app, Run from Signia is ideal for on-the-spot fittings. Incredibly comfortable thanks to its ergonomic shape, the convenient and easy-to-use Run delivers an outstanding hearing experience.

BTE			
Model	Run S	Run P	Run SP
Performance			
Matrix (2 ccm)	ThinTube 45/124 dB SPL Earhook 55/124 dB SPL	Earhook 70/132 dB SPL	Earhook 79/140 dB SPL
Battery			
Size	312	13	675
Power consumption	0.9 mA	1.3 mA	2.2 mA

Click ITE

Model		Run Click CIC	Run Click ITC
Performance			
Matrix (2 ccm)	50/113 dB SPL		55/116 dB SPL
Battery			
Size	10		312
Power consumption	1.1 mA		1.1 mA

Provide an outstanding hearing experience

- Key features
- A high level of amplification for severe to profound hearing loss
 - Ergonomic shape for a high level of wearing comfort
 - Simple fittings everywhere via the rocker switch

Basic

Fun

Equipped with Signia’s proven hearing aid technology, Fun covers all essential customer requirements. It offers superb sound quality for severe to profound hearing loss, excellent wearing comfort and is simply fitted by pushing the handy rocker switch.



Model	BTE	
	Fun P	Fun SP
		
Performance	Earhook 70/132 dB SPL	
Matrix (2 ccm)	Earhook 79/140 dBSPL	
Battery		
Size	13	675
Power consumption	1.2 mA	2.2 mA

Deliver a powerful listening experience

- Key features
- Simple fittings everywhere via the rocker switch
 - For mild to severe hearing loss
 - Easy handling

Basic

Fast

Fast BTE hearing aids offer a high level of performance and deliver just the right amount of amplification for customers with severe hearing loss. They are easily adjusted via the rocker switch for simplified fitting. Sleek and comfortable with a range of features, Signia’s Fast ensures an outstanding hearing experience suitable for every customer.



Model	BTE
	Fast P
	
Performance	
Matrix (2 ccm)	Earhook 62/130 dB SPL
Battery	
Size	13
Power consumption	1.1 mA

App, accessories and hearing aid compatibility

This overview shows you which apps and accessories work with our hearing aids. Depending on the solution your customers choose, they can even enhance their hearing experience and convenience with this well-matching portfolio.

Hearing aid	Signia app	CROS/ BiCROS	Streamers and transmitters	Remote Controls	Portable Charger	Inductive Charger II	Standard Charger	Dry&Clean Charger	Multi Charger
									
Pure Charge&Go IX	●	● ^N	●	●	●	-	●	●	●
Pure Charge&Go T IX	●	● ^N	●	●	●	-	●	●	●
Silk Charge&Go IX	●	● ^O	-	●	●	-	-	-	-
Styletto AX	●	● ^I	●	●	●	-	-	●	-
Insio Charge&Go AX	●	● ^J	●	●	-	-	●	-	-
Pure Charge&Go AX	●	● ^J	●	●	●	-	●	●	●
Pure Charge&Go T AX	●	● ^J	●	●	●	-	●	●	●
Pure 312 AX	●	● ^J	●	●	-	-	-	-	-
Active Pro	●	-	●	●	●	-	-	-	-
Motion Charge&Go X	●	● ^K	●	●	-	●	-	-	-
Motion Charge&Go P X	●	● ^K	●	●	-	-	●	●	●
Motion Charge&Go SP X	●	● ^K	●	●	-	-	●	●	●
Styletto X	●	-	●	●	●	-	-	-	-
Silk X	●	● ^L	-	-	-	-	-	-	-
Pure Charge&Go X	●	● ^K	●	●	-	●	-	-	-
Pure 312 X	●	● ^K	●	●	-	-	-	-	-
Motion 13 Nx	●	● ^C	●	●	-	-	-	-	-
Motion 13P Nx	●	● ^C	●	●	-	-	-	-	-
Insio Nx	●	● ^{G,H}	● ^F	● ^F	-	-	-	-	-
Intuis 3 RIC312	●	-	-	-	●	-	-	-	-
Intuis 4	●	● ^M	● ^M	● ^M	●	-	-	-	-
Intuis 3	●	-	-	-	●	-	-	-	-
Prompt	●	-	-	-	●	-	-	-	-
Run	●	-	-	-	-	-	-	-	-

C) Compatible with CROS Pure 312 Nx or CROS Pure Charge&Go Nx

D) Compatible with the Pure Charge&Go/Motion Charge&Go/Cellion charger

E) Compatible with CROS Silk Nx

F) Compatible with Bluetooth-enabled Insio Nx ITC/ITE

G) Insio Nx IIC and CIC are compatible with CROS Silk Nx

H) Bluetooth-enabled Insio Nx ITC and ITE are compatible with CROS Pure 312 Nx

I) Compatible with CROS Styletto AX

J) Compatible with CROS Pure Charge&Go AX and CROS Pure 312 AX

K) Compatible with CROS Pure Charge&Go X and CROS Pure 312 X

L) Compatible with CROS Silk X

M) only 4.1 and 4.2

N) Compatible with CROS Pure Charge&Go IX

O) Compatible with CROS Silk Charge&Go IX

Accessory overview

Streamers



Via StreamLine TV, TV sound can be fed directly into Bluetooth-equipped Signia hearing aids. Wearers can even adjust the TV volume individually via the Signia app.



The StreamLine Mic transforms Signia devices with Bluetooth into more than state-of-the-art hearing aids. They also become hands-free mobile headsets and wireless stereo headphones that stream highest quality audio into both ears from any Bluetooth device. Wearers can even adjust the streaming volume individually via the Signia app.

Remote Controls



Small enough to carry on a keyring, the miniPocket is the ideal solution for anyone who wants discreet, convenient hearing aid control without the need for a large remote control or smartphone.

CROS/BiCROS



For people with unaidable hearing loss on one side, Signia CROS devices receive sound on that side and transmit it to a hearing aid on the other side.

Charging Options



The Multi Charger is your ultimate one-stop charging solution. One charger compatible with a wide range of rechargeable hearing aids such as: Motion Charge&Go P X, Motion Charge&Go SP X, Pure Charge&Go AX, Pure Charge&Go T AX, Pure Charge&Go IX and Pure Charge&Go T IX for end users.



The chargers recharge hearing aids so your clients don't need to worry about handling batteries. The Inductive Charger II for Motion Charge&Go X and Pure Charge&Go X also dries hearing aids while charging to give a full day's ultra-durable hearing performance.



The Dry & Clean chargers for Pure Charge&Go IX, Styletto AX, Pure Charge&Go AX, Motion Charge&Go P X and Motion Charge&Go SP X offer a UV-C cleaning functionality while charging and drying as an additional benefit.



EarWear 3.0

Discreet. Secure. Comfortable.

EarWear 3.0 delivers superior comfort and fit. Its small, subtle design ensures a secure and barely noticeable presence. With easy-to-exchange Eartips and Sleeves 3.0, and compatibility with custom earmolds, it provides an optimal hearing experience for your clients


Your EarWear 3.0 advantages at a glance:

- High physical fit rate**
The combination of small receiver cross sections, soft eartips providing a narrow receiver front, and thoughtfully designed cable geometry in a range of lengths delivers a good fit for almost every ear.
- Wide range of acoustic coupling options**
Eartips and sleeves are available in a variety of styles, sizes, and venting options to ensure comfort and a good acoustic coupling.
- Good retention**
Cable geometry is designed to anchor the receiver reliability in the ear canal with no lateral migration, reducing the need for concha locks.
- Cosmetically appealing**
Clever cable shape and in-ear position ensure a close fit on the ear, reducing visibility.
- Wearer convenience**
Eartips and sleeves are easy for wearers to exchange when needed.
- Durable and hypoallergenic**
miniReceiver 3.0 is designed with a reinforced cable outlet for durability and a Titanium Shield for hypoallergenic wear.




EarWear 3.0

miniReceiver 3.0



- 1) push-/pull system – no tool required
- 2) Smart Receiver Detection prevents misuse
- 3) Unique curvature ensures high retention and sits close to head for cosmetic appeal
- 4) Five lengths 0-4
- 5) Reinforced cable outlet for high durability
- 6) Enhanced wax protection
- 7) Eartip 3.0 and Sleeve 3.0 compatible with miniReceiver 3.0 and ThinTube 3.0
- 8) WaxGuard 3.0 NanoCare
- 9) Interface allows easy Eartip 3.0 and Sleeve 3.0 exchange while offering high safety
- 10) Small interface enables comfortable fit in small ear canals
- 11) Hypoallergenic Titanium Shield



miniReceiver 3.0 available in four power levels: S, M, P, and HP


	S		M		P		HP	
	2ccm coupler	Ear simulator	2ccm coupler	Ear simulator	2ccm coupler	Ear simulator	2ccm coupler	Ear simulator
OSPL 90 (Peak)	110 dB SPL	120 dB SPL	119 dB SPL	129 dB SPL	122 dB SPL	131 dB SPL	131 dB SPL	138 dB SPL
138 dB SPL	46 dB	57 dB	60 dB	70 dB	65 dB	75 dB	75 dB	83 dB



EarWear 3.0

ThinTube 3.0

A choice of two diameters provides a cosmetically appealing solution for mild to severe hearing loss



Eartips 3.0 and Sleeves 3.0


This versatile new range is suitable for the vast majority of ear anatomies and hearing needs

Soft eartip and sleeve noses allow for deep insertion with high comfort – also into small ear canals

The secure and proven design makes them much easier for consumers to exchange.

The Power Sleeve 3.0 expands the acoustically highly reliable sleeve to also offer a fully closed option.


All Earwear 3.0 items can be used together.



EarWear 3.0

Earmolds 3.0

Compatible with miniReceiver 3.0 and ThinTube 3.0



Earmold 3.0 is available in a short and a long version.

The long Earmold 3.0 comes with:

- WaxGuard 3.0 NanoCare
- QuickGuard Wax Guard

Concha Lock Sleeve 3.0

Optional Concha Lock Sleeve 3.0 is available in three sizes



Concha Lock Sleeve 3.0 S – fits miniReceiver 3.0 S and both ThinTube 3.0 and ThinTube 3.0 P

Concha Lock Sleeve 3.0 M – fits miniReceiver 3.0 M

Concha Lock Sleeve 3.0 P – fits miniReceiver 3.0 P



EarWear 3.0

Tools and starter kits


Starter kits are available for each solution: Eartips 3.0, ThinTube 3.0 and miniReceiver 3.0.

Each kit includes specific items required for that type of EarWear 3.0. These include:


- Measurement gauge to ensure a perfect match for miniReceivers 3.0 and ThinTubes 3.0
- Matchbook containing 8 Wax Guard 3.0 NanoCare
- 10-pack of new one-size-fits-all Concha Locks 3.0
- Optional tool for exchanging eartips
- Tool to press out the receiver from a custom eartip

You also receive instruction manuals on how to use EarWear 3.0, covering everything from eartips and sleeves to attaching or removing ThinTubes 3.0 or miniReceivers 3.0 to exchanging earmolds or WaxGuards.


Starter kit Eartips 3.0 contents



Starter kit ThinTube 3.0 contents



Starter kit miniReceiver 3.0 contents



Removal Tool 3.0

EarWear 3.0

EarWear compatibility table

While newer Signia products use EarWear 3.0, some older products use EarWear 2.0

	ThinTube 3.0	miniReceiver 3.0	Eartips 3.0 & Sleeves 3.0	Earmold 3.0	ThinTube	Life-Tips	Life-Earmold	miniReceiver 2.0	Click Domes & Click Sleeves	Click Mold
Pure IX series	–	●	●	●	–	–	–	–	–	–
Silk IX series	–	–	● *	● **	–	–	–	–	–	–
Styletto AX series	–	●	●	●	–	–	–	–	–	–
Pure AX series	–	●	●	●	–	–	–	–	–	–
Motion X series	●	–	●	●	● ***	–	–	–	–	–
Pure X series	–	–	–	–	–	–	–	●	●	●
Styletto X series	–	–	–	–	–	–	–	●	●	●
Silk X series	–	–	–	–	–	–	–	–	●	–
Motion Nx series	–	–	–	–	●	●	●	–	–	–
Silk Nx series	–	–	–	–	–	–	–	–	●	–


* Silk Charge&Go IX compatible with Sleeve 3.0 only
** Custom earmolds for Silk IX series not interchangeable with standard Earmold 3.0
*** Motion SP X is not compatible with ThinTube



Meeting all your needs

The Signia hearing care professional support

Every day you try to meet and exceed your clients' hearing needs and wishes. In Signia you have a partner who fully supports you. Our broad range of hearing aid components, counseling and trial support as well as fitting tools enhances your workflow and provides you with all you need to find the best individual solution for your customers.



Signia Assistant

With its live artificial intelligence deep neural network, the Signia Assistant ensures that the end-user always gets the best possible solution for any given situation, always tailored to their specific needs and preferences.

It provides your clients with easy answers to all questions on handling, including showing them helpful videos in the app that are relevant to their specific device. In this way, people get the support they need during the entire onboarding process and you can spend the time at follow-up appointments more effectively.



TeleCare

Signia TeleCare enables you to support your clients remotely, so they get your help in real-life situations when they need it. You can hold virtual appointments via phone or video call and adjust the hearing aid settings remotely in the Connexx fitting software.

Reducing effort, increasing satisfaction

Support for counseling, fittings and trials



Connexx 9
fitting software



Unity™ 3
mobile fitting unit system

Signia offers you a broad range of counseling and fitting support – from programming software and tools to mobile fitting units to hearing aids for demonstration and advisory tools – up to comprehensive remote hearing care. This enhances your workflow and makes trials and fittings smoother and more pleasant for your customers.

Connexx with Integrated Xperience

Our Connexx fitting software puts you in the driver's seat with real-world insights to optimally support each individual patient.

First Fit in Connexx configures Signia IX hearing aids to allow the wearer to hear what is important to him or her. The First Fit settings have been carefully optimized to provide a good starting position for each wearer.

The Dynamic Soundscape Processing (DSP) control in Connexx provides a slider to adjust to your clients' individual preferences. This slider is now more powerful than ever thanks to our new Augmented Focus™ technology.

For you clients to fully optimize their personal sound in the situations important to them, they can use the Signia Assistant. The Assistant can fine-tune the DSP slider and show you if any changes were made from the initial settings. You can view details about each specific interaction: The situation your client was in at a particular date and time, the issue that required assistance, and the solution implemented by the Assistant.





All your clients' needs, all together


Signia app

Greater convenience

Simple and convenient, the Signia app empowers wearers with access to remote control adjustments, audio streaming, and the Signia Assistant. And because it includes TeleCare, it integrates seamlessly into your fitting process, making it easier for your customers to stay in touch with you when they are away from your premises. Suitable for use with all Signia hearing aids, its smart user interface only shows the user functions applicable to their own devices.



Compatible with Android and iOS.
The app is free of charge.



TeleCare

Stay in closer contact with your customers:

- Real-time contact via text message, voice call or video call
- Intuitive monitoring of your customers' daily satisfaction
- Set virtual appointments for the precise remote adjustment of the hearing aids

Remote Control

Personalize your customers' hearing experience further:

- Hearing program and volume control
- Sound Balance Adjustment
- 360° Spatial Configurator for individual adjustment of directionality
- Continuous monitoring of the hearing aids' connection and battery status

Connectivity & Streaming

- Easy management of audio streaming accessories (e.g. StreamLine TV, StreamLine Mic) as well as CROS devices
- Volume level control while streaming
- Adaptive streaming volume



TeleCare



Remote Control



Connectivity & Streaming



Signia Assistant


Signia Assistant

The smart new way to better address your clients' needs

Imagine being able to help your clients to Be Brilliant every single step of their hearing journey outside your office? Now you can – with the revolutionary Signia Assistant.


The world's first live deep neural network for hearing aids

Easy to use, the Assistant offers small sound adjustments when needed and how-to videos to answers handling questions. With its live artificial intelligence deep neural network, the Assistant gives your clients the best possible solution when issues arise in between appointments. All they need is their Signia IX or AX hearing aids and the Signia app.



See every detail with Connexx

With the Connexx fitting software, you'll find all the changes made by the Assistant, enabling you to understand each patient's key challenges. It's the smart solution for more efficient consultations and fittings - and more satisfied customers.



With the Signia Assistant you can:

- Give your clients the assurance of 24/7 support
- Increase your client satisfaction and conversion rate with the power to address hearing aid wearers' needs more accurately
- Increase the quality of follow up appointments through a new level of real-world usage insights
- Fine-tune with ease by seeing exactly how your clients choose to personalize their hearing experience

Color overview SLIM-RIC, RICs and BTEs

[illegible]

1) Styletto X 1/2 PL available in black or white only

Color overview Click ITEs

Standard colors for Click ITEs

Description	Type	Active Pro	Silk IX	Silk X	Intuis 3	Prompt	Run
 Blue/Mocha	left Click CIC		●	●			
 Red/Mocha	right Click CIC		●	●			
 Blue/Black	left Click CIC		●	●			
 Red/Black	right Click CIC		●	●			
 Mocha/Beige	Click CIC				●	●	●
 Tan/Beige	Click ITC				●	●	●
 Black	Click ITE	●					
 Black/Silver	Click ITE	●					
 Snow White/Rose Gold	Click ITE	●					

Further colors available on request.

Color overview ITEs



Color Flipper
Presenting RIC, BTE and ITE colors,
and product-specific overview of
available colors per product.

Item No. 03287-99C1-7600

Standard colors for custom-made ITEs

Description	Type	Insio C&G AX ITC/ITE	Insio Nx CIC / IIC	Insio Nx ITC/ITE	Intuis 3	Prompt
 Beige	ITC, ITE	● ¹	●	●	●	●
 Tan	ITC, ITE	● ¹	●	●	●	●
 Mocha	CIC, ITC, ITE	● ¹	●	●	●	●
 Brown	CIC, ITC, ITE	● ¹	●	●	●	●
 Dark Brown	IIC, ITC, ITE	● ¹	●	●	●	●
 Black	IIC, CIC, ITE	● ¹				

Further colors available on request.

(1) Available as ITC and ITE



Detailed feature overview SLIM-RIC and RICs

Type	Pure Charge&Go IX	Pure Charge&Go T IX	Styletto AX	Pure Charge&Go AX	Pure Charge&Go T AX	Pure 312 AX
Performance levels	7IX/5IX/3IX		7AX/5AX/3AX/2AX/1AX			
Picture						
Product details						
Matrix (2 ccm)	S 46/110 dB SPL M 60/119 dB SPL P 65/122 dB SPL HP 75/131 dB SPL	S 46/110 dB SPL M 60/119 dB SPL P 65/122 dB SPL HP 75/131 dB SPL	S 46/110 dB SPL M 60/119 dB SPL P 65/122 dB SPL	S 46/110 dB SPL M 60/119 dB SPL P 65/122 dB SPL HP 75/131 dB SPL	S 46/110 dB SPL M 60/119 dB SPL P 65/122 dB SPL HP 75/131 dB SPL	S 44/108dB SPL M 58/116 dB SPL P 62/119 dB SPL HP 74/128 dB SPL
Battery size	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	312
Power consumption						1.5-1.7 mA
Features						
Frequency channels/ Compression controls	7IX – 48/20 5IX – 36/18 3IX – 32/16 - -	7IX – 48/20 5IX – 36/18 3IX – 32/16 - -	7AX – 48/20 5AX – 32/16 3AX – 24/12 - -	7AX – 48/20 5AX – 32/16 3AX – 24/12 2AX – 16/8 1AX – 16/8	7AX – 48/20 5AX – 32/16 3AX – 24/12 2AX – 16/8 1AX – 16/8	7AX – 48/20 5AX – 32/16 3AX – 24/12 2AX – 16/8 1AX – 16/8
Hearing programs	6	6	6	6	6	6
e2e wireless 4.0	●	●	●	●	●	●
Direct streaming (via Mfi and ASHA)	●	●	●	●	●	●
TwinPhone	● ¹⁶	● ¹⁶	● ¹⁵	● ¹⁶	● ¹⁶	● ¹⁶
TeleCare	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁹
Tinnitus therapies (available on 7/5/3/2 PL)	●	●	●	●	●	●
CROS/BiCROS (available on 7/5/3/2/1 PL)	● ³⁰	● ³⁰	● ²⁹	● ²¹	● ²¹	● ²¹
Rocker switch	●	●	-	●	●	●
Volume control	● ²⁴	● ²⁴	-	● ²⁴	● ²⁴	● ²⁴
Rechargeable	●	●	●	●	●	-
T-Coil	-	●	-	-	●	-
Battery door on/off function	-	-	-	-	-	●
Audio input	● ²⁸	● ²⁸	● ²⁸	● ²⁸	● ²⁸	● ²⁸
Modifications						
Color conversion kit	○	○		○	○	○
App/Accessory compatibility						
Signia app	●	●	●	●	●	●
Smart Connect	-	-	-	-	-	-
StreamLine TV/StreamLine Mic	●	●	●	●	●	●
miniPocket	●	●	●	●	●	●
Connexx Smart Remote	-	-	-	-	-	-
ConnexxAir	-	-	-	-	-	-
Noahlink Wireless	●	●	●	●	●	●
Programming adapter/cable	-	-	-	-	-	-
Fitting set	miniReceiver 3.0	miniReceiver 3.0	miniReceiver 3.0	miniReceiver 3.0	miniReceiver 3.0	miniReceiver 3.0

(15) To be activated via Signia app or miniPocket
(16) To be activated via rocker switch, Signia app or miniPocket
(19) Full live remote tuning via Connexx
(21) With either CROS Pure 312 AX or CROS Pure Charge&Go AX

(24) Programmable via Connexx
(28) Via direct streaming
(29) Available on 7/5/3 PL with CROS Styletto AX
(30) Available on 7/5/3 PL with CROS Pure Charge&Go IX

● Standard
○ Optional



Type	Styletto X	Pure Charge&Go X	Pure 312 X	Intuis 3 RIC312
Performance levels	7X/5X/3X/2X/1X			Essential
Picture				
Product details				
Matrix (2 ccm)	S 45/108 dB SPL M 60/119 dB SPL P 70/124 dB SPL -	S 45/108 dB SPL M 60/119 dB SPL P 70/124 dB SPL HP 75/130 dB SPL	S 45/108 dB SPL M 60/119 dB SPL P 70/124 dB SPL HP 75/130 dB SPL	S 45/108 dB SPL M 60/119 dB SPL P 70/124 dB SPL HP 75/130 dB SPL
Battery size	Lithium-ion power cell	Lithium-ion power cell	312	312
Power consumption			1.7 – 1.9 mA	0.9 – 1.1 mA
Features				
Frequency channels/ Compression controls	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	12/6
Hearing programs	6	6	6	4
e2e wireless	● ¹³	● ¹³	● ¹³	-
Direct streaming	●	●	●	-
TwinPhone	● ¹⁵	● ¹⁶	● ¹⁶	-
TeleCare	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁸
Tinnitus therapies (available on 7/5/3/2 PL)	●	●	●	-
CROS/BiCROS (available on 7/5/3 PL)	-	● ²¹	● ²¹	-
Rocker switch	-	●	●	● ²³
Volume control	-	● ²⁴	● ²⁴	● ²³
Rechargeable	●	●	-	-
T-Coil	-	-	○	●
Battery door on/off function	-	-	●	●
Audio input	● ²⁸	● ²⁸	● ²⁸	-
Modifications				
Color conversion kit	-	○	○	○
App/Accessory compatibility				
Signia app	●	●	●	●
Smart Connect	-	-	-	-
StreamLine TV/StreamLine Mic	●	●	●	-
miniPocket	●	●	●	●
Connexx Smart Remote	-	-	-	
ConnexxAir	-	-	-	-
Noahlink Wireless	●	●	●	-
Programming adapter/cable	-	-	-	size 312
Fitting set	miniReceiver 2.0	miniReceiver 2.0	miniReceiver 2.0	miniReceiver 2.0

(13) Ultra HD e2e with binaural audio signal processing
(15) To be activated via Signia app or miniPocket
(16) To be activated via rocker switch, Signia app or miniPocket
(17) To be activated via push button, Signia app or miniPocket
(18) Basic remote tuning via 4-band-equalizer

(19) Full live remote tuning via Connexx
(21) With either CROS Pure 312 Nx or CROS Pure Charge&Go Nx
(23) Exchangeable controls – programmable via Connexx
(24) Programmable via Connexx
(28) Via direct streaming
(29) Push button

● Standard
○ Optional



Detailed feature overview BTEs

Type	Motion Charge&Go X		Motion Charge&Go P X		Motion Charge&Go SPX		Motion 13 Nx		Motion 13P Nx		
Performance levels	7X/5X/3X/2X/1X						7Nx/5Nx/3Nx/2Nx/1Nx				
Picture											
Product details											
Matrix (2 ccm)	Earhook ThinTube/ LifeTube	60/133 dB SPL 3.0 60/125 dB SPL -	77/135 dB SPL 3.0 66/130 dB SPL -	82/140 dB SPL - -	63/134 dB SPL 56/126 dB SPL -	77/135 dB SPL 65/130 dB SPL -					
Battery size		Lithium-ion power cell	3.0 66/130 dB	Lithium-ion power cell	13	13					
Power consumption			-		1.4 – 2.0 mA	1.7 – 2.0 mA					
Features											
Frequency channels/ Compression controls		7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7X – 34/16 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8					
Hearing programs		6	6	6	6	6					
e2e wireless		● ¹³	● ¹³	● ¹³	● ¹³	● ¹³					
Direct streaming		●	●	●	●	●					
TwinPhone (available on PL7/5/3)		● ^{2, 16}	● ^{2, 16}	● ^{2, 16}	● ¹⁶	● ¹⁶					
TeleCare		● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁹					
Tinnitus therapies (available on 7/5/3/2 PL)		●	●	●	●	●					
CROS/BiCROS		● ²	● ²	● ²	● ^{2, 21}	● ^{2, 21}					
Push button		-	-	-	-	-					
Volume control		-	-	-	● ²⁴	● ²⁴					
Rechargeable		●	●	●	-	-					
T-Coil		-	●	●	○	○					
Battery door on/off function		-	-	-	●	●					
Audio input		● ²⁸	● ²⁸	● ²⁸	● ²⁸	● ²⁸					
Modifications											
Color conversion kit		○	○	○	○						
App/Accessory compatibility									○		
Signia app		●	●	●	●	●					
Smart Connect		-	-	-	-	-					
StreamLine TV/StreamLine Mic		●	●	●	●	●					
miniPocket		●	●	●	●	●					
easyPocket		-	-	-	-	-					
ConnexxAir		-	-	-	-	-					
Noahlink Wireless		●	●	●	●	●					
Programming adapter/cable		-	-	-	size 13	size 13					
Fitting set		Thin Tube 3.0/Earhook	Thin Tube 3.0/Earhook	Earhook	ThinTube/Earhook	ThinTube/Earhook					

(2) Not available with PL 2/1
(13) Ultra HD e2e
(16) To be activated via rocker switch, Signia app or miniPocket
(18) Basic remote tuning via 4-band-equalizer
(19) Full live remote tuning via Connexx
(24) Programmable via Connexx
(28) Via direct streaming
(33) With battery compartment modification

Intuis M 4	Intuis P 4	Intuis SP 4	Intuis 3 S	Intuis 3 M	Intuis 3 P	Intuis 3 SP
Essential						
60/133 dB SPL 3.0 61/125 dB SPL -	77/135 dB SPL 3.0 66/129 dB SPL -	82/140 dB SPL - -	55/124 dB SPL - 45/124 dB SPL	60/130 dB SPL - 53/125 dB SPL	70/134 dB SPL - -	80/140 dB SPL - -
13	13	6,75	312	13	13	675
1.6 mA	1.9 mA	2.2 mA	0.9 mA	1.0 mA	1.4 mA	2.4 mA
4.2 – 16/8 4.1 – 16/8 4.0 – 12/6	4.2 – 16/8 4.1 – 16/8 4.0 – 12/6	4.2 – 16/8 4.1 – 16/8 4.0 – 12/6	12/6	12/6	12/6	12/6
4	4	4	4	4	4	4
●	●	●	-	-	-	-
● ⁴³	● ⁴³	● ⁴³	-	-	-	-
-	-	-	-	-	-	-
●	●	●	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸
● ⁴⁴	● ⁴⁴	● ⁴⁴	-	-	-	-
● ⁴³	● ⁴³	● ⁴³	-	-	-	-
-	-	-	● ²⁴	-	●	●
●	●	●	-	● ²⁴	● ²⁴	● ²⁴
-	-	-	-	-	-	-
●	●	●	-	●	●	●
●	●	●	●	●	●	●
-	-	-	-	○ ³³	○ ³³	○ ³³
●	●	●	○	○	○	-
●	●	●	●	●	●	●
-	-	-	-	-	-	-
● ⁴³	● ⁴³	● ⁴³	-	-	-	-
●	●	●	●	●	●	●
-	-	-	-	-	-	-
-	-	-	-	-	-	-
●	●	●	-	-	-	-
-	-	-	size 312	size 13	Programming shoe/ CS 44 cable	Programming shoe/ CS 44 cable
ThinTube/Earhook	ThinTube/Earhook	ThinTube/Earhook	LifeTube	LifeTube	-	-

(43) only 4.1 and 4.2
(44) only 4.2
● Standard
○ Optional



Detailed feature overview BTEs

Type	Prompt S	Prompt P	Prompt SP	Run S
Performance level	Basic			
Picture				
Product details				
Matrix (2 ccm)	ThinTube/ LifeTube Earhook	-	-	-
	45/124 dB SPL - 55/124 dB SPL	- 70/134 dB SPL	- 80/140 dB SPL	45/124 dB SPL - 55/124 dB SPL
Battery size	312	13	675	312
Power consumption	0.9 mA	1.4 mA	2.4 mA	0.9 mA
Features				
Frequency channels/ Compression controls	8/4	8/4	8/4	8
Hearing programs	4	4	4	6 preconfigured sound profiles
e2e wireless	-	-	-	-
Direct streaming	-	-	-	-
Push button	●	●	●	-
Volume control	-	● ⁽²⁴⁾	● ⁽²⁴⁾	●
T-Coil	-	●	●	-
Battery door on/off function	●	●	●	●
Audio input	-	○ ⁽³³⁾	○ ⁽³³⁾	-
App/Accessory compatibility				
Signia app	●	●	●	●
Smart Connect	-	-	-	-
miniPocket	●	●	●	●
Programming adapter/cable	size 312	Programming shoe/ CS 44 cable	Programming shoe/ CS 44 cable	MobileFitting app
Fitting set	LifeTube	Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling

(24) Programmable via Connexx
(33) With battery compartment modification

● Standard
○ Optional



Run P	Run SP	Fun P	Fun SP	Fast P
Basic				
				
- - 70/132 dB SPL	- - 79/140 dB SPL	- - 70/132 dB SPL	- - 79/140 dB SPL	- - 62/130 dB SPL
13	675	13	675	13
1.3 mA	2.2 mA	1.2 mA	2.2 mA	1.1 mA
8	8	6	6	4
6 preconfigured sound profiles	6 preconfigured sound profiles	4 preconfigured sound profiles	4 preconfigured sound profiles	3 preconfigured sound profiles
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
●	●	●	●	●
-	-	-	-	-
●	●	●	●	●
-	-	-	-	-
●	●	-	-	-
-	-	-	-	-
-	-	-	-	-
MobileFitting app	MobileFitting app	Fitting via rocker switch	Fitting via rocker switch	Fitting via rocker switch
Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling	Simple acoustic coupling



Detailed feature overview Click ITEs and ITEs

Type	Silk Charge&Go IX Click CIC	Insio Charge&Go AX	Active Pro	Active	Silk X Click CIC
Performance levels	7IX/5IX/3IX	7AX/5AX/3AX/2AX/1AX			7X/5X/3X/2X/1X
Picture					
Product details					
Matrix (2 ccm)	50/114 dB SPL	S 50/113 dB SPL ¹ M 55/118 dB SPL ¹ P 65/124 dB SPL ¹	42/118 dB SPL	42/118 dB SPL	50/114 dB SPL
Battery size	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	10
Power consumption					1.2 mA
Features					
Frequency channels/ Compression controls	7IX – 48/20 5IX – 36/18 3IX – 32/16	7AX – 48/20 5AX – 32/16 3AX – 24/12 2AX – 16/8 1AX – 16/8	48/20	16/8	7X – 48/20 5X – 32/16 3X – 24/12 2X – 16/8 1X – 16/8
Hearing programs	6	6	6	4	6
e2e wireless	● ¹³	● ⁴¹	● ¹³	● ¹³	● ¹³
Direct streaming	-	●	●	●	-
TwinPhone (available on PL7/5/3)	● ¹⁵	● ¹⁵	● ¹⁵	-	● ¹⁵
TeleCare	● ¹⁸	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁸
Tinnitus therapies (not available with PL 1)	●	●	●	-	●
CROS/BiCROS	●	● ⁴⁰	-	-	●
Push button	-	●	-	-	-
Volume control	-	-	-	-	-
T-Coil	-	-	-	-	-
Battery door on/off function	-	-	-	-	●
Audio input	-	● ³⁹	● ²⁸	● ²⁸	-
Optivent	-		-	-	-
App/Accessory compatibility					
Signia app	●	●	●	●	●
Smart Connect	-	-	-	-	-
StreamLine TV/StreamLine Mic	-	●	●	●	-
miniPocket	●	●	●	●	●
easyPocket	-	-	-	-	-
ConnexxAir	-	-	-	-	●
Noahlink Wireless	-	●	●	●	-
Programming adapter/cable	Flex connector				Flex connector
Fitting set	Sleeve 3.0 / EarTip 3.0 / Earmold 3.0	-	Click Sleeve 2.0 / Click Dome	Click Sleeve 2.0 / Click Dome	Click Sleeve 2.0 / Click Dome

(1) Max. values
(13) Ultra HD e2e with binaural audio signal processing
(15) To be activated via Signia app or miniPocket
(18) Basic remote tuning via 4-band-equalizer
(19) Full live remote tuning via Connexx
(21) With CROS Pure 312 Nx
(25) Standard on Bluetooth-enabled ITC/ITE using
Signia app/miniPocket, optional on other ITC/ITE
(28) Via direct streaming
(29) With Smart Connect

(34) Standard on Bluetooth-enabled ITC/ITE,
optional on other ITC/ITE
(35) To be activated via Signia app, miniPocket
or optional push button
(37) Not available for IIC
(38) Only if e2e wireless is available
(39) Through Bluetooth
(40) With either CROS Pure 312 AX or CROS Pure Charge&Go AX
(41) e2e wireless 4.0
(42) not for IIC

● Standard
○ Optional



Insio Nx IIC/CIC	Insio Nx ITC	Insio Nx ITE	Intuis 3 Click CIC	Intuis 3 Click ITC	Intuis 3 CIC	Intuis 3 ITC	Intuis 3 ITE
7Nx/5Nx/3Nx/2Nx/1Nx			Essential		Essential		
							
S 50/113 dB SPL ¹ M 55/118 dB SPL ⁴² P 65/124 dB SPL ⁴²	S 50/113 dB SPL ¹ M 55/118 dB SPL ¹ P 65/124 dB SPL ¹	M 55/118 dB SPL ¹ P 65/124 dB SPL ¹	50/113 dB SPL	55/116 dB SPL	50/113 dB SPL 55/118 dB SPL 65/124 dB SPL ¹	50/113 dB SPL 55/118 dB SPL 65/124 dB SPL ¹	55/118 dB SPL 65/124 dB SPL ¹
10	312	312	10	312	10	312	13
1.3 mA ¹	1.3 mA ¹	1.3 mA ¹	1.1 mA	1.1 mA	1.1 mA ¹	1.1 mA ¹	1.1 mA ¹
7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	7Nx – 48/20 5Nx – 32/16 3Nx – 24/12 2Nx – 16/8 1Nx – 16/8	12/6	12/6	12/6	12/6	12/6
6	6	6	4	4	4	4	4
○ ¹³	○ ¹³	○ ¹³	-	-	-	-	-
-	● ³⁴	● ³⁴	-	-	-	-	-
● ³⁵	● ³⁵	● ³⁵	-	-	-	-	-
● ¹⁸	● ^{19, 34}	● ^{19, 34}	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸
●	●	●	-	-	-	-	-
●	● ²¹	● ²¹	-	-	-	-	-
○ ³⁷	○ ³⁴	○ ³⁴	-	●	○	○ ³⁴	○ ³⁴
-	-	-	-	-	-	-	-
-	-	-	-	-	-	○	○
●	●	●	●	●	●	●	●
-	○ ²⁸	○ ²⁸	-	-	-	-	-
○	○	○	-	-	○	○	○
●	●	●	●	●	●	●	●
-	-	-	-	-	-	-	-
-	● ³⁴	● ³⁴	-	-	-	-	-
●	●	●	●	●	●	●	●
-	-	-	-	-	-	-	-
● ³⁸	-	-	-	-	-	-	-
-	●	●	-	-	-	-	-
Flex connector ³⁷	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector
-	Click Sleeve 2.0 / Click Dome	-	Click Sleeve 2.0 / Click Dome	Click Sleeve 2.0 / Click Dome	-	-	-



Detailed feature overview Click ITEs and ITEs

Type	Prompt Click CIC	Prompt Click ITC	Prompt CIC	Prompt ITC	Prompt ITE
Performance levels	Basic				
Picture					
Product details					
Matrix (2 ccm)	50/113 dB SPL	55/116 dB SPL	50/113 dB SPL 55/118 dB SPL 65/124 dB SPL ¹	50/113 dB SPL 55/118 dB SPL 65/124 dB SPL ¹	55/118 dB SPL ¹ 65/124 dB SPL ¹
Battery size	10	312	10	312/10	13/312
Power consumption	1.1 mA	1.1 mA	1.1 mA ¹	1.1 mA ¹	1.1 mA ¹
Features					
Frequency channels/ Compression controls	8/4	8/4	8/4	8/4	8/4
Hearing programs	4	4	4	4	4
Push button	-	●	○	○	○
Volume control	-	-	○ ²⁵	○ ²⁵	○ ²⁵
T-Coil	-	-	-	○	○
Battery door on/off function	●	●	●	●	●
Optivent	-	-	○	○	○
App/Accessory compatibility					
Signia app	●	●	●	●	●
miniPocket	●	●	●	●	●
Programming adapter/cable	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector
Fitting set	Click Sleeve 2.0 / Click Dome	Click Sleeve 2.0 / Click Dome	-	-	-

(1) Max. values
(25) Standard on Bluetooth-enabled ITC/ITE using Signia app/miniPocket,
optional on other ITC/ITE

● Standard
○ Optional

Run Click CIC	Run Click ITC
Basic	
50/113 dB SPL	55/116 dB SPL
10	312
1.1 mA	1.1 mA
8	8
6 preconfigured sound profiles	6 preconfigured sound profiles
-	-
-	-
-	-
●	●
-	-
●	●
-	-
MobileFitting app	MobileFitting app
Click Sleeve 2.0 / Click Dome	-