

6 REASONS

to choose Virtuo Vivo™







Scanning in real color



Sterilizable sleeves



Small & lightweight



Open STL files



Streamlined workflows

dental wings

New! Virtuo Vivo™ Remarkably compact and powder-less.



Resembling a standard dental tool, the small and lightweight (7.5 oz / 213 g) wand ensures unparalleled patient comfort and allows the clinician to assume a natural position relative to the patient when scanning.



Design

Mill

Scan

SMARTER DENTISTRY

Streamlined workflows that take scans even further.

Chairside Design & Production

Transfer files for chairside design with DWOS® Chairside or with other 3rd party software.



Real color scanning allows for a precise and user-friendly design of restorations.



Mill with the Roland milling machines

Surgical Planning & Production

Use intraoral impressions for surgical planning and integrated drill guide design with coDiagnostiX®.



Send files directly to coDiagnostiX or export the STL files to a 3rd party surgical planning software.



Orthodontics

Connect directly to ClearCorrect™ clear aligner provider for orthodontic treatment.



Clear aligners are effective, removable, and patients are increasingly opting for this solution.



Outsourcing to Lab

With open STL files, send digital impressions directly to a lab via DWOS Connect or other export method for design and manufacturing.







Voor meer informatie: Ditare B.V. 076- 587 6969 www.ditare.info info@ditare.info